

# SleepRight<sup>®</sup> NO-BOIL<sup>®</sup> DENTAL GUARD

## TEETH GRINDING PROTECTION & TMJD PREVENTION

No boiling required for fitting • Most durable bitepad material available • 17-years of proven safety


### HERE'S HOW IT WORKS!

The SleepRight Rx's adjustable bitepads, made of a hybrid resilient material, fit comfortably between your teeth and are designed to re-establish your natural freeway space. While the SleepRight Rx protects teeth from clenching and grinding, we strongly encourage our patients to concentrate on, "lips together, teeth apart" to initiate the muscle relaxation learning process.

The SleepRight Rx does not require boiling water. Detailed pictures are included with easy to follow care and fitting instructions. With SleepRight Rx's patented adjustable technology, you have four bitepad positions, making it easy to customize the size to fit most people. The bitepads adjust forward for a smaller fit and backward for a larger fit, as well as swivel to match the angle of all bites.

The SleepRight Rx utilizes natural tissue retention between your lips, cheek, and tongue to be held in place. The flexible, heat sensitive band is positioned between your lower lip, cheek and gums and uses your body's temperature to quickly adapt to each individual's anatomy.

### SELF-ADJUST TECHNOLOGY



### Dura-Comfort<sup>®</sup>

- Flexible stability wings for extra safety
- Frontal band is anatomically designed for more comfort



### Slim-Comfort<sup>®</sup>

- Ideal for small mouths and sensitive gums
- Slim design for comfort



### 30-DAY MONEY BACK GUARANTEE

This guarantee allows you to try the SleepRight Rx with your patients and evaluate it for yourself. If for any reason you are not satisfied, simply return it to us within 30-days from the original purchase date for a refund or credit towards your next purchase. Please include proof of purchase.

### 1-YEAR REPLACEMENT WARRANTY

This warranty begins when you give the product to the patient. If the patient grinds through it within the first 12-months, you may give the patient a SleepRight Rx from your inventory, send the used one back to us for either a replacement or credit towards your next purchase. This policy also applies if the patient would like to exchange the Dura-Comfort<sup>®</sup> for a Slim-Comfort<sup>®</sup> or vice versa, or if the product simply doesn't work for the patient. Patient may receive a replacement up to 2 times within the year. Please include proof of purchase.

INSURANCE COVERAGE: D9940 - Occlusal appliance for bruxism D7880 - Occlusal orthotic device

1.888.792.0865 • [www.sleepright.com/professional](http://www.sleepright.com/professional)